

Latest generation digital radiography solutions with  
**Integrated, Dynamic, Portable, Wired and Wireless detectors**



Réalisation : A.S. éditions - 01 77 55 05 42 - www.beckitions.com - STEPHANIX se réserve le droit de modifier les caractéristiques de ses produits en fonction de l'évolution perpétuelle des technologies. Photographies non contractuelles. Crédits Illustration : l'atelier 06 80 10 22 83 - STEPHANIX autorise le droit de modifier les produits. Features regarding the evolution of the technologies without notice. Non-binding document. D114

# OMNISCOP SERIES



**Head Office**  
 Rue Jean Moulin - ZI du Bayon  
 42150 La Ricamarie - FRANCE  
 Tel : + 33 4 77 47 81 60  
 contact@stephanix.com

**stephanix**  
 RADIOLOGICAL SOLUTIONS

**www.stephanix.com**

**OMNISCOP SERIES**  
**Export department**  
 15, rue du Saule Trapu  
 91882 Massy cedex - FRANCE  
 Tel : + 33 1 60 13 56 25  
 Fax : + 33 4 77 37 55 19

**stephanix**  
 RADIOLOGICAL SOLUTIONS

SEE THOROUGHLY, ENSURE SAFETY

# OMNISCOP SERIES

## SYSTEM OVERVIEW

**OMNISCOP SERIES** is positively impact your outcomes and enhances your productivity. This system is equipped with a high frequency generator from 5 to 15 kW, fixed or rotated anode and an image intensifier from 9" to 12" image intensifier that handles a complete range of procedures : surgery, traumatology, orthopedics, vascular\* and cardiac\*.

The small footprint of the C-arm facilitates patient positioning and working position.

The high-resolution CCD camera performs high image quality at the lower dose for an optimum definition according the procedures. The post-processing highlights tiny details and produces the

clear and accurate images clinical demands. As data management is important, **OMNISCOP SERIES** can be connected to hospital network thanks to DICOM interface.

**OMNISCOP SERIES** is an ideal system for ER, OR, ICU, surgical imaging or interventional setting.

\*depending on the model



## VERSATILITY SOLUTION

**OMNISCOP SERIES** is designed to give you more freedom. The small footprint ensures easy positioning of the unit around the operating table.

All movements of the C-arm are balanced for a fast and accurate positioning.

Easy to drive handles for side to side movement and smooth transport.



## EASE OF USE

**OMNISCOP SERIES** is equipped with a comprehensive and intuitive user interface and open structure workstation.

The touch screen console allows obtaining an accurate image. A short list of APR is preset, the operator selects the protocol adjusts the imaging parameters for a warranty of the best image quality and lower dose levels.

Ergonomically shaped, the workstation is designed for an optimal visibility in operating room. The image processing ensures a quick access for patient data and post-processing is easy.

